

**Amendments to the Specification**

Please replace the paragraph beginning at line 13 of page 13 as follows:

The core member of the different sized and shaped US-ICD will all be the same size and shape. That way, during an implantation ~~procedures~~ procedure, multiple sized US-ICDs can be available for implantation, each one without a core member. Once the implantation procedure is being performed, then the correct sized US-ICD can be selected and the core member can be inserted into the US-ICD and then programmed as described above. Another advantage of this configuration is when the battery within the core member needs replacing it can be done without removing the entire US-ICD.